

what's your vial style?

get accurate TOC results with vials for all your applications

Whether you're testing for compliance, cleaning validation, or difficult-to-recover compounds, Sievers* TOC vials ensure top performance and minimize risks. Accurate, efficient, & confident – that should be your style.

- Certified, traceable, and quaranteed
- Manufactured by the global leader in TOC
- Backed by our OOS investigation support



Sievers Certified (<10 ppb TOC) Vials

for reliable, guaranteed quality

- · The right choice for regulatory compliance testing
- Obtain accurate TOC results the first time and every time
- Pre-cleaned and certified (<10 ppb TOC)
- Backed by our OOS investigation support



Dual Use Conductivity & TOC (DUCT) Vials

for simultaneous compendial testing

- Suitable for simultaneous Stage 1 conductivity and TOC compliance testing
- Using a single vial and automation saves time, eliminates sample handling issues, and enhances data integrity
- Specialty coated glass and specialty septa, no ionic leaching, certified (<10 ppb TOC)
- · Backed by our OOS investigation support



Pre-Acidified Vials

for sticky proteins and peptides

- Improve protein recoveries for cleaning validation applications
- Prevent adherence of proteins and peptides to vial surfaces
- Available in options for swab samples (pre-filled with acidified water) or for final rinse samples
- Backed by our OOS investigation support



Certified Blue Vials

for cleanrooms

- Used in controlled environments
- Paper free packaging, produced in an ISO 7 Cleanroom
- Certified Blue packaging available for Sievers Certified, DUCT, and Pre-Acidified Vials
- Backed by our OOS investigation support



Trademark of SUEZ: may be registered in one or more countries